

Dr. Schumacher GmbH · Am Roggenfeld 3 · 34323 Malsfeld · GERMANY

14 April 2025

2025-343

Rinsing after ULTRASOL OXY® WIPES

To Whom It May Concern:

With reference to your request, whether it is necessary to rinse medical device surfaces after using our ULTRASOL OXY® WIPES, we are able to provide the following information:

After the drying time, surfaces/medical devices that have been cleaned and disinfected with ULTRASOL OXY® WIPES can be used immediately. Our wipes leave no toxic residues. This makes the product user friendly and it can also be used in the food sector.

In addition, these products leave no toxic or environmentally harmful residues on surfaces. After drying, the active ingredients decompose into oxygen, water and a very light acetic acid (0.1 %), so that skin irritations according to EN ISO 10993-23 (2021); part 23 (skin irritation human epidermis) and EN ISO 10993-10 (2023); part 10 (skin sensitization) can be excluded. The residues after drying were tested with reference to the impregnating solution according to EN ISO 10993-5 (2023); Part 5 (In vitro cytotoxicity tests) and classified as non-cytotoxic.

Further information regarding our products can be found on our homepage under:

<https://www.schumacher-online.com/en/downloads/data-sheets/disinfection-hygiene/>

If you have any further questions, please do not hesitate to contact us.

Kind regards

Dr. Schumacher GmbH



Christian Frank
Director of Customer Science



Beatrix Witteler
Key Account Manager International Sales